

KEN BEDPAN WASHER DISINFECTOR 738

Intelligent Disinfection

CLEAN, LEAN & GREEN

- ✓ Total disinfection
- ✓ Easy service access
- ✓ Complies with 15883

The KEN BWD 738 has been developed for hospitals, special wards and nursing homes. It is useful for washing and disinfection of bedpans, urine bottles, washing bowls, suction glasses and similar utensils.

With its compact design, it offers a large capacity. Furthermore, the machine is environmentally optimised with a low consumption of energy, water and detergents, which will result in low operational costs.

KEN BWD 738 is supplied as a put-through model with two doors, which offers a secure separation between the clean and unclean side.

At the unclean side the door can be operated manually or automatically by means of a foot pedal, and on the clean side the door is opened manually.

The machine is easy to operate by touch panel with icons showing all relevant information of the program needed by the operators.

The machine is functional with all types of bedpans, urinal bottles, commodes, urine glasses and a special utensils basket is available for washing kidney bowls and other small items.

All data from the BWD 738 are fully exchangeable with common computer systems.



HEALTHCARE



All KEN machines are in conformity with present EU directives, and all the hospital machines have been constructed in accordance with EN ISO 15883.

Advantages:

Environment

- Low energy costs
- Low noise level
- Low chemical usage

User-friendliness

- Easy to use icons
- Quick-start - automatic closing
- Door suitable as slopsink - bedpan emptied automatically

Hygiene safety

- EN 62366 tested
- EN ISO 15883-1-3-5
- Disinfection guarantee
- Separation between clean and unclean side
- Ao value programmable from 60 - 3000

| Description | Values |
|-------------------------------|---|
| Outer dimensions | H: 1745 x W: 760 x D: 560 mm |
| Dimensions washing chamber | H: 500 x W: 485 x D: 370 mm |
| Weight | 141 kg |
| Integrated combi tank | 9 L |
| Built in detergent cabinet | Room for 2 x 5 L cans |
| Capacity, alternatively- | 1 bedpan + 3 urine bottles OR 3 bedpans (2 with urine) OR 3 washing bowls OR 1 utensil basket OR 1 suction bottle OR 1 bucket |
| Bedpan program time | Appr. 6 min |
| Sealing class | IP 21 (not rinse-proof) |
| Water connections | Cold: 3/4" G (poss. 1/2" pipe), press.: 200-800kPa, temp.: 5-30°C, hardness: 0-3°dH Hot: 3/4" G (poss. 1/2" pipe), press.: 200-800kPa, temp.: 5-70°C, hardness: 0-3°dH |
| Sound level | Appr. 60 dB (A) |
| Water consumption per phase | See technical manual |
| Power supply options | 3N x 200 - 440V AC 50/60Hz +PE+(N) (Single phase available) |
| Heating element in combi tank | 6 kW / 9 kW |
| Programs available | 5 standard + 1 optional |
| Drain: S or P water trap | Inner ø90 mm, outer ø110 mm |

Contact a KEN partner

KEN HYGIENE SYSTEMS A/S
Bøgebjergvej 60 · DK-5672 Broby
T +45 6263 1091 · F +45 6263 1607
ken@ken.dk · www.ken.dk

Clean, Lean & Green

